



ISSUE NO	DATE	REVISION / ECO NUMBER	APP
1	1-16-84		JIC
2	4-16-84	ADDED ALARM OUT	JIC

TELEDYNE  
ANALYTICAL INSTRUMENTS  
JUN 1 2010  
DOCUMENT CONTROL  
COPY

- NOTES:**
- 1- ALLOW CLEARANCE IN FRONT OF ANALYZER FOR OPENING OF HOUSING DOOR (13 1/2" OR 343 MILLIMETERS)
  - 2- ALLOW CLEARANCE ON SIDES OF UNIT FOR TUBING & ELECTRICAL HOOK-UP; CLEARANCE NOT REQUIRED ON BACK OF SYSTEM.
  - 3- ALL DIMENSIONS IN INCHES ± 1/16" OR MILLIMETERS ± 1.0

ITEM	PART NO.	REQ	DESCRIPTION
TOLERANCE UNLESS OTHERWISE SPECIFIED		THIS DRAWING IS THE PROPERTY OF TELEDYNE ANALYTICAL INSTRUMENTS AND CONTAINS CONFIDENTIAL INFORMATION. IT IS NOT TO BE COPIED, REPRODUCED OR USED WITHOUT WRITTEN PERMISSION.	
FRAC. 1-64	TELEDYNE ANALYTICAL INSTRUMENTS CITY OF INDUSTRY, CALIFORNIA 91743		SCALE QUARTER DATE 1.13.84
DEC. 005			MAT'L N
ANGULAR . . .			FINISH N
MDL 328 R PERCENT OXYGEN ANALYZER OUTLINE DIAGRAM			DR. ADZANO ENGR APP
			B-33976