

# Fermenter Off Gas

## Biotechnology for Specialty Chemicals

### Purpose

Falcon GC analysis for real-time fermentation monitoring.

### Measured Components

- Chemicals in fermenter off-gas
- Measurement in ppm range

### Cycle Time

< 60 seconds

### Detector(S)

FID

### Reference Methods

N/A

\*Contact Teledyne For Specifications On This Application.



**TELEDYNE**  
**ANALYTICAL INSTRUMENTS**  
Everywhereyoulook™

16830 Chestnut Street, City of Industry, CA 91748 USA  
Phone: 626-934-1500/888-789-8168 Fax: 626-934-1651  
EMAIL: Ask\_tai@teledyne.com www.teledyne-ai.com

***Contact the TAI Sales Team to ask about our other instruments.***