

Food Grade Fatty Acids

Purpose

Improve batch process throughput with fast, automated, food grade, fatty acids end point analysis using the Teledyne-Falcon GC.

Measured Components

- C3 – C50 food grade fatty acids SIMDIS and discrete hydrocarbon analysis
- Measurement in ppm to %

Cycle Time

<12 minutes

Detector(S)


Dual FID


Reference Methods

N/A

*Contact Teledyne for specifications on this application.

Chromatograms





16830 Chestnut Street, City of Industry, CA 91748 USA
Phone: 626-934-1500/888-789-8168 Fax: 626-934-1651
EMAIL: Ask_tai@teledyne.com www.teledyne-ai.com

Contact the TAI Sales Team to ask about our other instruments.