

Sulfur Compounds in Liquid Phase Hydrocarbon Samples

Purpose

Lab analysis of discrete hydrocarbons in condensate samples at well heads and fractionators with Falcon GC.

Measured Components

- Hydrogen sulfide, methyl mercaptan, dimethyl sulfide, ethyl mercaptan, carbonyl sulfide
- Measurement range ppm

Cycle Time

360 seconds (6 minutes)

Detector(S)

FPD

Reference Methods

N/A

*Contact Teledyne For Specifications On This Application.



16830 Chestnut Street, City of Industry, CA 91748 USA
Phone: 626-934-1500/888-789-8168 Fax: 626-934-1651
EMAIL: Ask_tai@teledyne.com www.teledyne-ai.com

Contact the TAI Sales Team to ask about our other instruments.