



Ultrafast Gas Chromatography in Transportable and OnSite Applications

John Crandall, presenter, Ned Roques, Matt Holliday





Abstract



Ultrafast Gas Chromatography in Transportable and Onsite Applications

John Crandall, presenter, Ned Roques, Matt Holliday – Falcon Analytical Brian Rohrback – Infometrix

J.C. Reyes, Bill Winniford – The Dow Chemical Company

Transportable and onsite advanced analytical chemistry especially gas chromatography or gas chromatography/mass spectrometry usually leads to giant motor home style platforms. These systems use either enormous shore power cables or electrical generators taking hours to set up and to bring the analytical systems up to stable operating condition. Often, the need is rapid and possibly even stealthy operational capability. With social media, vehicles can easily divert their route to avoid inspection stations. Examples include roadside detection of markers used at the ppm level in untaxed fuel or adulteration of fuels with potentially engine damaging materials to avoid tax.

Implementation of systems in the UK and Republic of Ireland will be discussed as well as prospects for expansion throughout the European Union and the rest of the world. Prospects for deterring the adulteration of fuel and even counterfeit adult beverages "in the back of a Jeep" will be discussed.




CALIDUStm Intellectual **Property**




Fast Micro Gas Chromatograph System US Patent 8414832



Trans-Configurable Modular Chromatographic Assembly US Patent 8336366









How does CALIDUS Ultrafast GC work? **Electronics** 4444 Hypothetical Sample Mixture сссс a = gasoline b = jet Vent **Detector** c = diesel **Signal** abcbcacab bcacababc cababcbca abcbcacab bcacababo cababcbca abcbcacab **Stainless Steel Capillary Gas Chromatography Separation Column: the heart of the system** Thin Film or Micropacked Technology **Carrier Gas** + Source (~50 psi)

The Calidus Modular GC System



- Sample Processing Module
 - Split/splitless injector
 - Septum purge
 - Inlet glass liner
 - Column switching
- Column Modules
 - 2 16 meters for 32 meters total
 - 180 530 micron ID
 - Various film thickness
 - Even micropacked available
- Detectors
 - Flame Ionization
 - Thermal Conductivity
 - Flame Photometric
 - Dielectric Barrier Discharge with Helium Ionization, Electron Capture & Photo Ionization modes





Swappable Choice of Detectors


Choice of Columns

Flowing Gas or Liquid Sample GC System




Septum Purge


Modules Combined into Instrument Models




101 – single column 2m – 16m, single detector




201 - single column 4m - 32m, single detector in series




301 – dual column 2m – 16m, dual detector in parallel



CS- dual column 2m – 16m, dual detector with column switching valve, up to 32 meters total



GS – single column 2m – 16m, with large sample loop, single detector with sample pump



GS-T – single column 2m – 16m, with preconcentration trap, single detector with sample pump

New



Fuel Laundering Examples (UK/Ireland)

The problem: removal of "red dye" from agricultural and marine fuels saves road use tax... at the expense of "the people's roads."

Republic of Ireland (March 2013)



Laundering by bleaching agent

UK (March 2011)



Laundering by acid treatment

Probably NOT environmentally friendly!



Dow Fuel Marking: Usage, History, and Value

□ Usage

 Fuel markers are used by government agencies and fuel marketers to differentiate fuel sources

□ Value

- Global fuel dyes and markers market estimated at \$342 MM by 2019
- In 2013 UK and Ireland lost an estimated 1.7 billion of tax revenue as a result of illegal fuel laundering activities

☐ Current Dow Marking Products

- AUTOMATE™ and MORTRACE™ colored fuel dyes
- SPECTRACE™ digital fuel markers
- ACCUTRACE™ launder resistant forensic fuel markers

□ ACCUTRACE™ S10: Key Benefits


- ☐ most <u>resistant</u> to all known criminal removal techniques
- □ highly cost effective in deterring fuel laundering, which funds criminal activity and results in environmentally hazardous waste dumping
- ☐ introduces no known ozone-depleting additives as it is made of same elements as fuel itself
- ☐ <u>invisible</u> to the eye and difficult for criminals to discover but <u>readily detectable</u> by law enforcement using highly accurate and mobile testing methods



Heartcut 2D GC Method

□ 2D GC with heartcutting capability using a Deans switch allows separation of components that are unresolved on the primary column

- ☐ avoids false positives due to co-elution of other components present in the fuel matrices.
- ☐ heartcut 2D GC identified as method of choice for detection and quantification.







CALIDUS™ CS, for roadside analysis of ACCUTRACE™



- ☐ highly transportable
 - ~11 kilograms
 - < 300 Watts
 - ~ 17 x 8.5 x 11 inches
- **□**economic
 - lowest capital
 - lowest consumables
- ☐ highest Durability
 - field proven
 - industry standard FIDs
 - temperatures well below spec limit on the CS valve














CALIDUS™ CS Instrument Layout




☐ Columns and Detector

- o relatively non-polar column on Channel 1
- o relatively polar column on Channel 2
- o rugged, reliable, industry standard FID detectors
 - FID 1 used to determine heartcut timing and normalization of results
 - FID 2 used to measure ACCUTRACE™ S10

☐ CS Valve for Heartcut

- high temperature rotary valve with pneumatic actuation
- special temperature control employed on valve body
- temperature well below the specification limit for the valve






Resistively heated column modules

Calibration Standard on Separation Channel 1




Expanded time scale of Initial Separation Channel




Material Heartcut into MXT-50 including ACCUTRACE™ S10


Calibration Standard on Heartcut Channel 2






Linearity / LOD, LOQ





Triplicate injections at $49 \mu g/L$ (overlaid with unmarked diesel; red trace)

ACCUTRACE™S10 in diesel



LOD ppb (μg/L)	19
LOQ ppb (μg/L)	63
LOQ equivalent to dosage rate	0.63%

- □ ACCUTRACE™ S10 quantifiable down to 63 ppb (µg/L).
- ☐ Setting a trigger at 125 ppb and a tolerance value insures no false negatives and no false positives.



Implementation in the Field (UK/Ireland)



Fuel Marker Demonstrates Success in United Kingdom & Ireland

Governments Report Progress in Fight Against Fuel Laundering





"Since the launch of the new marker we have seen a significant reduction in the detection of laundering plants and dumping of waste material associated with laundering activities."
- Pat Curtis¹⁰

National Oils Lead Investigator, HMRC

- in the first six months, HMRC made 26,261 tests with 1,299 detections made.
- 85 cases identified where ACCUTRACE[™] S10 was present with low levels of 'old style' markers (indication of laundering attempts)
- increased fuel sales of legitimate diesel by 2.1% (1Q 2014) and 5.9% (1Q 2015).
- fuel duty on legitimate diesel increased by £310 million in the first 3 months of use as illegal fuel launderers ceased or decreased their operations.

European Union completing their Health, Safety and Environmental study now

- 2 African nations contemplating projects in 2018
- S. Korea and Brazil have programs developing

Award Winning Team



- ☐ A robust 2D-GC heartcut portable detection method was developed to distinguish marker from fuel matrix.
- □ ACCUTRACE™ S10 currently used in marking tax rebated fuels in UK and Ireland with successful implementation of CALIDUS™ ultrafast GC in the field.
- **□** Evaluation report by HMRC and IGEES showed that the use of Dow's ACCUTRACE™ S10 Marker has led to:
 - reduction in laundering plants discovered; a decline in illegal use of subsidized fuels; a drop in waste products associated with fuel laundering.
 - negligible selling of laundered fuel and the problem is close to being eliminated.
- **□** Dow's ACCUTRACE[™] fuel marking & detection technology recognition:
 - 2016 R&D 100 Awards
 - 2016 IChemE Global Awards
 - Edison Awards
 - Presidents E Awards
 - WV Export Awards












Pipeline Interface Detection... Interference Free N₂ @ 35 Seconds













The Result: 170 Second Resolution for Determining Flare Start & Stop Timing






OnSite N₂ Analysis




Before

Grab Sample Analysis Profile




After
Near Real-time Analysis
Profile



Far less wasted Product 1 & 2 and Nitrogen

Fully Transportable System for Designer Diesel,






EU Diesel, 10% Base Oil Added and 5% Base Oil 3% Veg Oil

Same as Previous, Y-Axis Expanded





In the Hands of Fuel Inspectors



Program Launch

Dedicated user interface for analytical and data security

chromperfect

Fuel Sentry


Set up the acquisition hardware: 1. Turn on the hydrogen generator. Press Open when 100%. 2. Switch on the Calidus GC. 3. Press Next to initialize the Calidus GC.

Step by Step Instructions

Reduces needs for printed SOP

On Screen Results

Dependable results on screen in about 3 minutes



Automatic Daily Calibrations

Assures instrument is working within specifications every day

Data Input Fields

Quit Next

Analysis of RT mix standard is complete.


Press Next to make another run. Press Quit to shut down the GC.

*** PASS ***

Field notes integrated with data file


Blank Run






RT Standard





RGO, Lot#2





RGO, Lot#2



Reference Gas Oil - QC Report

RGO Data File: 0214171020_[R]GO.0149.CDF Blank File: 0214171008_[B]lank.0147.CDF


Reference Gas Oil QC Detail - Type: Ref Oil Lot #2


Percent Off	Actual BP Temp(C)	Temp Window	Actual BP Temp(C)	Result
IBP	115.0	7.00	115.16	-PASSED-
10	176.1	4.44	175.64	-PASSED-
20	223.9	5.00	224.73	-PASSED-
30	259.4	4.78	259.96	-PASSED-
40	288.9	4.28	289.42	-PASSED-
50	312.2	4.28	312.58	-PASSED-
60	331.7	4.28	332.48	-PASSED-
70	353.9	4.28	354.27	-PASSED-
80	377.8	4.28	378.49	-PASSED-
90	406.7	4.28	407.67	-PASSED-
FBP	475.0	11.78	471.94	-PASSED-





"Designer Diesel"























"Designer Diesel"























"Designer Diesel"




















Adult Beverages, Food Ingredients & Contamination


Food & Beverage relating to product quality, authentication and contamination




Beverage

- Ethanol and fusel oils in wine and distilled spirits
- Validation that ethylene glycol is NOT present (Italy)
- Soon in Hungary, forensic authentication for wine & liquor





Ethanol & Fusel Oils in Distilled Spirits
 Jose Cuervo, Grey Goose, Johnnie Walker
 Black, & Crown Royal



Questions?





Thank you for your interest & attention!

