

Model 2000A / 2000AEU

Orders should include the part number, model number, and instrument serial number. Call Teledyne Analytical Instruments at +1-888-789-8168 (USA Only) or +1-626-934-1500 to place an order for spare parts. Contact your regional sales representative for instrumentation questions; see [http://ai.www-
uat.teledyne.com/en-us/who-we-are/worldwide-representatives](http://ai.www-uat.teledyne.com/en-us/who-we-are/worldwide-representatives).

Recommended Parts		
QTY	Part No.	Description
1	C-65507A	Back panel board
1	C-62371A	Front panel board
1	C-65098	Preamplifier board
1	C-73870D	Main computer board
1	CP1802	Shielded cable clamp
50	CP1799	Solder cup contact for CP1798 (pins)
1	CP1798	PLU, 50 pin D-Sub connector
1*	F9	Fuse, 1A, 250V, 3AG, Slowblow (US)
2*	F1275	Fuse, 1A, 250V, 5 x 20mm, T-Slowblow (EU)
1	A-33748	Thermistor assembly
1	B-66119	34 pin connector ribbon cable assembly, front to main PCB
1	B-67065	14 pin connector ribbon cable assembly LCD to main PCB
1	A-31157	Standard cell assembly
1	A-35013	SS cell block with gold filaments
* Order one type only, US or EU as appropriate		



16830 Chestnut Street, City of Industry, CA 91748 USA
Phone: 626-934-1500/888-789-8168 Fax: 626-934-1651
EMAIL: Ask_tai@teledyne.com www.teledyne-ai.com

Find your local sales representative: [http://ai.www-
uat.teledyne.com/en-us/who-we-are/worldwide-representatives](http://ai.www-uat.teledyne.com/en-us/who-we-are/worldwide-representatives)