

Model 2020

Orders should include the part number, model number, and instrument serial number. Call Teledyne Analytical Instruments at +1-888-789-8168 (USA Only) or +1-626-934-1500 to place an order for spare parts. Contact your regional sales representative for instrumentation questions; see [http://ai.www-
uat.teledyne.com/en-us/who-we-are/worldwide-representatives](http://ai.www-uat.teledyne.com/en-us/who-we-are/worldwide-representatives).

Recommended Parts		
QTY	Part No.	Description
1	C73870-D	Main PC Board assembly
1	B-67674	Integrally mounted SS cal/zero control valves
1	F-448	Flowmeter (reference)
1	F-449	Flowmeter (sample)
1	D-67472	Back panel / Interface board (includes temp control circuit)
1	A-31157	Standard cell assembly
1	A-35013	SS cell block with gold filaments
1	A-33748	Thermistor assembly
1	H-158	Heater assembly
1	A-73957-F	Eprom



16830 Chestnut Street, City of Industry, CA 91748 USA
Phone: 626-934-1500/888-789-8168 Fax: 626-934-1651
EMAIL: Ask_tai@teledyne.com www.teledyne-ai.com

Find your local sales representative: [http://ai.www-
uat.teledyne.com/en-us/who-we-are/worldwide-representatives](http://ai.www-uat.teledyne.com/en-us/who-we-are/worldwide-representatives)